

Scalable Production of hiPSC-Derived Cardiomyocytes in Stirred-Tank Bioreactors

Philipp Nold^{1*}, Karin Langenberg², Viola Borgdorff², Marta Duńska², Pieter van Loenen², Farbod Famili², Yang Jiang², Arie Reijerkerk^{2*} ¹Eppendorf AG Bioprocess Center, Juelich, Germany; ²Ncardia BV, Leiden, The Netherlands

*Contacts: nold.p@eppendorf.com; arie.reijerkerk@ncardia.com

Abstract

Stem cell-based technologies lay the basis for pioneering approaches in drug discovery, drug safety testing, and regenerative medicine. The routine application of human induced pluripotent stem cells (hiPSCs) and functional cells derived from hiPSCs in high-throughput applications for drug discovery will require the constant supply of high cell numbers at a consistent, high quality. Two-dimensional systems like T-flasks are widely used for the cultivation of stem cells; however, they are limited in control and scalability. Stirred-tank bioreactors have emerged as promising culture system for large-scale expansion of stem cells and their differentiation into the desired cell type. Ncardia[®] is a stem cell drug discovery and development company developing highly predictive human cellular (disease) assay systems for safety and efficacy testing. The company's goal is to accelerate and improve drug candidate selection, reduce drug development costs and ultimately increase drug discovery and development efficiency. Largescale manufacturing of hiPSC-based cardiomyocytes at high quality is integral for its cardiovascular services. Here we demonstrate the suitability of the Eppendorf DASbox[®] Mini Bioreactor system for the expansion of hiPSCs and subsequent production of Ncardia's Pluricyte[®] Cardiomyocytes. We show that using the DASbox system, hiPSCs were successfully expanded as cell aggregates in a highly reproducible manner. During the production of three different batches, aggregates were of homogenous size and the cells retained key iPSC markers including Sox2, Oct3/4 and Nanog during the expansion phase. More importantly, using the DASbox, a fully controlled production process scalable for large-scale manufacturing of hiPSC-derived cardiomyocytes could be established. Quantification of cardiomyocyte markers (cTNT and MLC2v) revealed the robust and efficient generation of multiple and pure ventricular-like cardiomyocyte batches. Analysis of the electrophysiological properties by multielectrode arrays (MEA) demonstrated the mature phenotype of the cells. The expected response to cardioactive reference compounds further confirmed the functionality of bioreactor-derived cardiomyocytes. Altogether, the Eppendorf DASbox Mini Bioreactor System allows smooth integration of Pluricyte Cardiomyocytes into high-throughput cardiovascular screening solutions to enhance the assessment of cardiac safety and efficacy of drug candidates.

Stem cell technology for drug discovery

The data presented were obtained by scientists at Ncardia. They tested, whether the DASbox Mini Bioreactor System was suitable for hiPSC expansion and differentiation in a reproducible manner.





Figure. 1. Ncardia aims to use human stem cell technology to get better medicines to patients faster by improving the drug discovery and development process.

Bioprocess equipment

Figure. 2: DASbox Mini Bioreactor

Marker expression in DASbox

Figure 3. Characterization of iPSC aggregates and





System

The DASbox Mini Bioreactor System was used as culture system for the expansion of hiPSCs as cell aggregates and their differentiation into cardiomyocytes.
Parallel control of up to 24 bioreactors
Allows monitoring and control of critical process parameters like dissolved oxygen, pH, temperature
For use with conventional glass bioreactors or

 parameters like dissolved oxygen, pH, temperature
 For use with conventional glass bioreactors or BioBLU[®] Single-Use Vessels
 Rigid-wall, stirred-tank design simplifies culture scale-up

differentiated cardiomyocytes in DASbox Mini Bioreactors

show reproducible stem cell and cardiac marker expression.
A: Phase contrast imaging of iPSC aggregates
B: FACS profiles of iPSC marker expression at the start of pre-cardiac specification phase

C: Quantification of iPSC marker FACS data

D: Percentage of cTNT in cardiomyocytes batches generated with DASbox

Characterization of iPSC-derived cardiomyocytes



Figure 5. Expected pharmacological responses of bioreactor-derived cardiomyocytes to CiPA refer-

ence compounds

A: Isoproterenol (β-adrenergic receptor agonist) reduces beat period (enhances beat rate). Mexilitine blocks the rapid inward sodium current as well as the hERG potasium channels leading to a reduction in sodium spike amplitude (**B**,**C**) and a prolongation of the field potential duration (**D**,**E**) respectively (**C**,**E** show averaged waveforms). Nifedipine (calcium channel blocker) causes a shortening of the field potential duration (**F**). hERG potassium channel blockers dofetilide (**G**) and E4031 (**H**) prolong the field potential duration (**G**,**H**) and cause flattening of the repolarization peak (**I**,**J** show average waveforms). Arrhythmia's were observed for both dofetilide (10 nM, **K**) and E4031 (30 nM, **L**). Corresponding DMSO concentrations were used as control.







0.3-

0.2

Figure 4. Phenotypic and functional characterization of three batches of iPSC-

derived, cryopreserved cardiomyocytes produced in bioreactors A: Percentage of cTNT+ and

B: cTNT+/MLC2v+ co-expression in three cryopreserved batches of bioreactor-derived Pluricyte Cardiomyocytes.

C: Phase contrast image of thawed cells, 8 days after culture in Pluricyte Cardiomyocyte Medium.

D: Functionality of the cardiomyocytes (sodium spike amplitude, beat period, beat period irregularity, field potential duration, number of active electrodes per well; >300 µV sodium spike and the number of FPD electrodes; with a detectable T-wave) was assessed in Maestro[™] MEA system at different timepoints after plating.

High-throughput screening

All	1	2 3 4	5	6	7 8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	DR
A	MM	MMrmmmm	MMM	nmul	MINIM	un nu	MMMM	MMM	MALLAN	www.	MANN	www	MMM	MAMAN	www	MMMM	IMAAMA	non	mm	nn	-	3 n
в	MMM	mmmmmmm	NNNN	Inner	MMM	www	NUMBER		NNNN	222	-1-1	M	·/	MMMM		MUMUM	Internation	~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~			-	10 n
С	mm	mmmmm	MMM	mm	MMM	WW AND		,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	MMM	-h-		-r	-1-		NANGANANA	ananan	MANAMAN	and and a				30 n
D	0.1%	DMSOMMAN	humm	wwwww	MMMM	MMM mm			MANAN	1	1-	r-1										<u>100 n</u>
E	MM	MAAAAAAAA	aqueanantag	www.weiweiw	ANA	umu WMI.		AMMMA	MIM	-		1-				-						300 n
F	WWW	MAMMMAN				IIIII		human	1000000	1-	-1-		-	_								<u>1μ</u>
G	-l-l	MANNANANA					ANATTANA			J.J.	M	~~	-					_				3 µ
н	0.3 µ	MBayK 8644	DRC	com	pound	1 D	RC com	pour	nd1/2/1	DRC	com	pour	nd 3	DRC	com	pour	nd 4	DRO	C con	npour	nd 5	10 µ
1	summents	and the second second second	warmer.			MM	mmm	MMM	mmm	WWWW	www	MAA	M	MMA	NNW	MMM	MMM	MMM	MMM	mmm	www	
L							MMM	mm	MMM	MM	m	m	MM	MM	MMM	man	n	NNN	MMM	mm	mm	
К		water water to an and the second	w			NN	WWW	MMM	MMM	m	w	m	NN	MMM	MM	MW	www	MMM	hrmm	mm	www	
L	0.3 µ	MNitrendipine					MMMM	MANA	MM	mm	W	ww	nn	MM	-n-n	harland	mom	MANA	mm	manan	www	
М	mmm		1	-			-mm	MM	~~~	m	m	m	M	- Ma	-l-	N	~~	NUVAU	mm	MMM	MM	
Ν	MININ	MMMMMMMMMMMMMMMMMMMMMMMMMMMMMMMMMMMMMMM	Į	_		~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~	MM	~~~ v	mar	m	NN	NN	SV		r.			NUUN	man	www	www	
0		NNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNN					-vv	n	~~~~	MM	n-v-	s	MMM					mm	MAN	mm	MM	
P	ШµМ	Isoprotereno	DRC	com	pound	6 DI	RC com	pour	nd 7	DRC	com	pour	nd 8	DRC	com	npour	1d 9	No	treat	ment	MAN	

Figure 6. High throughput compound screening in

FLIPR calcium transient assay using bioreactor-derived

cardiomyocytes.

Bioreactor-derived cardiomyocytes were cultured for 10 days in 384-well plate format. On day 10 cells were treated with dose response curve (DRC) of various cardioactive compounds. 0.1 % DMSO was used as negative control and BayK 8644, Nitrendipine and Isoproterenol were used as positive controls. Effect of compounds 30 minutes post treatment were assessed in FLIPR calcium transient assay. The expected effects of compounds on calcium signal were observed in various concentrations as follows: Isoproterenol increased peak frequency, BayK 8644 increased peak amplitude and peak width, Nitrendipine decreased peak amplitude.

Conclusion and perspectives

- > DASbox Mini Bioreactor System was suitable for expansion of hiPSCs and differentiation into hiPSC-derived cardiomyocytes:
- > DASbox-derived cardiomyocytes showed excellent baseline electrophysiological function, revealed expected phenotypic responses to a set of reference compounds, and
- were suitable for high-throughput compound screening in a calcium transient assay.
- > Perspective: Further implement bioreactor-derived Pluricyte Cardiomyocytes into high-throughput cardiovascular screening solutions to enhance the assessment of cardiac safety and efficacy of drug candidates.

Acknowledgements

 This project has received funding from the European Union's Horizon 2020
 Research and Innovation Programme under grant agreement No. 726513.

