

# Optimize your cell culture and fermentation workflows with reliable, customizable solutions

## Cedex Bio and Cedex Bio HT Bioprocess Analyzers<sup>†</sup>

Count on high performance and reliability

- Broad dynamic range with onboard dilution
- Low, 200 µl sample volume
- Expandable assay menu with barcoded reagents
- Low maintenance



## MycotoOL PCR Mycoplasma Detection Kit<sup>†</sup>


Release your product in < 5 hours

- **FDA accepted** for biopharmaceutical release testing
- Sensitive to < 1 CFU/ml for all validated species
- Detection of more than 140 Mollicutes species
- No false positives, due to carryover prevention

## In vitro Glycoengineering Portfolio

Increase your glycosylation homogeneity

- High-quality enzymes and activated sugars
- Scalable model amenable to GMP
- Cost-effective solutions for discovery through commercialization



To learn how these products can improve your workflow, visit Roche Booth 309.

<sup>†</sup> For use in quality control/manufacturing process only.